

Thoughtfully designed for maximum convenience

23.8" bezel-less screen
for an immersive experience



Auto wake-up of residual
power reminder



iConsole—intelligent and
clinical exam-specific panel
layout with E-ink keys



A minimum height of 1m
Can be folded up and
taken on-the-go
Easily fits into an MPV



Robust flight case
Easy packaging



Dual-wing floating arm for
flexible adjustments

15.6" full-HD touch screen
with intuitive interaction



Full-space floating control
panel with an exceptionally
wide range of adjustments



26dB super-silent design as
quiet as a bedroom
Up to 4H uninterrupted
scanning without power cables



Modular design for easy
and fast maintenance

Resona I9

Elite Edition

Diagnostic Ultrasound System

Innovation, in every facet





Innovation

The **Resona I9** provides an entirely new experience, thanks to the innovations it features both inside and out. At its core is the revolutionary ZST+ platform, which delivers unparalleled image quality, infinite imaging solutions in dedicated applications, comprehensive multi-parametric quantification, and remarkable intelligence for improved efficiency.


Prepare to be greeted by an ultrasound system design that is truly one-of-a-kind. Take advantage of its latest innovative features to minimize your fatigue during scans.



Unlock infinite imaging with the Elite edition

Explore more details with UMA

Ultra Micro Angiography (UMA) is an innovative technology that surpasses the limitations of traditional Doppler imaging. By enhancing flow sensitivity and resolution, UMA revolutionizes diagnostic accuracy by offering expanded visibility into blood flows down to minute vessels, thereby boosting confidence in the diagnosis process.



Elevated imaging performance and expanded clinical adaptability



Convex



Linear



Phased array




Endocavity




Others

Superb confidence with unrivalled clarity


ABD



Liver HiFR STE




UWN+ CEUS




Fatty liver lab - USAT


Breast & thyroid



Smart Breast




Smart Thyroid




Breast mass STE


Cardiovascular



Carotid artery stenosis




RIMT




TT QA


Others



Bowel image



Smart Bladder



HR flow in metastatic lymph node